

Fissativo per Color Restore

Glossy transparent fixative, to be applied for finishing the treatment with Color Restore cod 45CR. Suitable for both natural and synthetic leather.
Shake well before use.

CHARACTERISTICS

- Physical state: liquid
- Base: water
- Color: transparent

APPLICATION

The product can be applied

- with a sponge like the one supplied;
- with a brush;
- with a spray gun, using EGO Carbonio 190 HALO Ø 0.5-0.7mm nozzle.

For each single application, we recommend two coats spaced 3 hours apart: the first light, the second more full-bodied.

Wait at least 8 hours for complete drying.

The effectiveness of our products is based on practical experiences and research work carried out in our laboratories; nevertheless we accept no liability for work carried out following our instructions being clear that the final result depends in all cases on a series of unforeseeable factors.

* For any information about product codes or packs, please see our catalogue, our price list or contact us.